

CADD™-Solis VIP Pump Guide

Medication/Rate: _____

If your medicine is stored in the refrigerator, take it out and let it warm to room temperature. **Do not heat the medicine in a microwave oven, direct sunlight, hot water, or any other heat source. Heating may damage the medicine and make it harmful.**

Follow these instructions with every infusion:

1. Prepare a clean work area.
2. Wash your hands well using a liquid soap. Be sure to wash all parts of your hands including the top and bottom, between fingers, and under fingernails. Use a paper towel to dry your hands. Use an alcohol-based hand sanitizer if no soap and water are available.
3. Supplies needed:
 - Medicine:
(CADD™ cassette with medicine CADD™ tubing set with medicine in an IV bag)
 - Pump
 - “Patient Information Guide” booklet
 - 4 AA Batteries
 - Alcohol pads
 - Saline Flushes (if instructed by the nurse or pharmacist)
 - Heparin Flush (if instructed by the nurse or pharmacist)
4. Follow the “Patient Information Guide” sections for starting your infusion:
 - a. Installing Batteries
 - b. Turning the Pump On and Off
 - c. Removing a Cassette
 - d. Attaching a New Cassette
 - e. Priming the Tubing
 - f. Resetting the Reservoir Volume
 - g. Starting the Pump
 - h. Stopping the Pump

For questions or concerns, please call Regional Home+ Home Infusion at (605) 755-1150.